



Alloys Catalogue





CONTENT

High gold bonding alloys

- 2 V-Gnathos® PF
- 2 V-Gnathos® Supra
- 2 V-Gnathos® Plus
- 3 V-Supragold
- 3 V 92
- 3 V-Classic

Medium gold bonding alloys

- 4 V-Deltaloy
- 4 V-Delta® Special
- 4 V-Delta® SF
- 5 V-Delta® 450

Palladium-based bonding alloy

- 5 Cerapall® 6
- 5 Cerapall® 2
- 6 Ceradelta® 2
- 6 Ceradelta®

Non precious alloys

- 6 Metalloy CC

Universal alloys

- 7 Aurofluid® Plus
- 7 Pontor® LFC
- 7 Pontor® 4 CF

High gold casting alloys

- 8 Aurofluid® 2 PF
- 8 Pontor® MPF
- 8 Aurofluid® 3

Medium gold casting alloys

- 9 Pontor® 2
- 9 Solaro 3
- 9 Yellow Special

Silver-palladium casting alloys

- 10 Pagalinor® 2
- 10 Pagalin® 2

Implant alloys

- 11 Implantalloy 1
- 11 Implantalloy 2
- 11 Implantalloy 3

Procedure

- 12 Casting
- 13 Surface preparing
- 14 Soldering

Solders for bonding alloys

- 16 VS 1 C
- 16 VS 1 Plus
- 16 VS 1 GPF
- 17 VS 1 A
- 17 Vacu PF
- 17 Vacu 2

Solders for universal alloys

- 18 Standard LFC 1

Solders for casting alloys

- 19 Vacu-PF
- 19 Standard S 4 PF
- 19 Standard S 2
- 20 Standard S 3
- 20 Standard S 31
- 20 Pagasol S 1

Other solders

- 21 Revetor II
- 21 Solchrom Pal

Solders for non precious alloys

- 22 Meta CC S1
- 22 Meta CC S2

Precious alloy wires

- 23 Ancrofil
- 23 Novostil

Gilding paste

- 24 Aurofilm® 2000
- 24 Aurofilm® NP

Flame chart

- 26 Hints for success

Oxidation table

- 28 Oxidation colour



V-GNATHOS® PF

Yellow, palladium-free, silver-free and copper-free, hard

V-Gnathos PF is the ideal alloy for esthetic single crowns. With its optimal balance between colour, hardness and melting range. It is easy to work with. It provides an esthetic yellow base for the application of porcelains.


CTE 14.6 / 14.9 10⁻⁶K⁻¹

Recommended solders

before bonding: VS 1 GPF / VS 1 Plus
after bonding: Vacu PF

Ordering information

200004784 25 g casting ingots
400003101 5 laser wires ø 0.3 mm



V-GNATHOS® SUPRA

Yellow, palladium-free, silver-free and copper-free, extra hard

V-Gnathos Supra is compatible with almost all high fusing ceramics, it combines a superior esthetics with an improved polishing property. V-Gnathos Supra is a guarantee for optimal patient safety.


CTE 14.3 / 14.5 10⁻⁶K⁻¹

Recommended solders

before bonding: VS 1 Plus
after bonding: Vacu PF

Ordering information

200004855 25 g casting ingots
400003446 5 laser wires ø 0.3 mm



V-GNATHOS® PLUS

Yellow, palladium-free, silver-free and copper-free, extra hard

The leading high gold bonding alloy. V-Gnathos Plus combines high hardness and high esthetic results for single crowns and short span bridges. It is unsurpassed in its biocompatibility.


CTE 14.3 / 14.6 10⁻⁶K⁻¹

Recommended solders

before bonding: VS 1 Plus
after bonding: Vacu PF

Ordering information

200004811 25 g casting ingots
400003102 5 laser wires ø 0.3 mm





V-SURPAGOLD


Light yellow, silver-free and copper-free, extra hard

V-Supragold is the modern alloy for long span bridges.
CTE 14.3 / 4.6 10⁻⁶K⁻¹



Recommended solders
before bonding: VS 1 GFP
after bonding: Vacu PF

Ordering information
200004825 25 g casting ingots
400003159 5 laser wires ø 0.3 mm



V 92


Light yellow, silver-free and copper-free, extra hard

Its melting range makes it easy to cast. V 92 is suitable for long span bridges, and has only one non-precious component. V 92 is an excellent milling alloy.
CTE 14.0 / 14.2 10⁻⁶K⁻¹



Recommended solders
before bonding: VS 1 A / VS 1 GFP
after bonding: Vacu PF

Ordering information
200004777 25 g casting ingots
400003110 5 laser wires ø 0.3 mm



V-CLASSIC


White, extra hard

V-Classic is the alloy of choice for use with attachments and precision milling. It is easy to solder and easy to polish. V-Classic has excellent thermal stability, high hardness and sag resistance.
CTE 14.1 / 14.4 10⁻⁶K⁻¹



Recommended solders
before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information
200004770 25 g casting ingots
400003106 5 laser wires ø 0.3 mm





V-DELTA[®]LOY

White, copper-free, extra hard

V-Deltaloy has a light oxidation colour and is easy to solder. Its wide range of indications makes V-Deltaloy a multi purpose alloy.


CTE 14.3 / 14.5 10⁻⁶K⁻¹

Recommended solders

before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information

200004807 25 g casting ingots
400003123 5 laser wires ø 0.3 mm



V-DELTA[®] SPECIAL

White, extra hard

V-Delta Special is an economical alloy due to a reduced content of palladium.


CTE 14.5 / 14.7 10⁻⁶K⁻¹

Recommended solders

before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information

200004832 25 g casting ingots
400003165 5 laser wires ø 0.3 mm



V-DELTA[®] SF

White, silver-free and copper-free, extra hard

V-Delta SF is versatile with excellent working, soldering and polishing properties. Its high hardness and high thermal stability make V-Delta SF a superb alloy for long span bridges.


CTE 13.8 / 14.1 10⁻⁶K⁻¹

Recommended solders

before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information

200004747 25 g casting ingots
400003111 5 laser wires ø 0.3 mm





V-DELTA® 450


White, extra hard

V-Delta 450 delivers excellent value and all-round performance.
CTE 14.0 / 14.2 10⁻⁶K⁻¹



Recommended solders
before bonding: VS 1 A
after bonding: Vacu PF

Ordering information
200004815 25 g casting ingots
400003161 5 laser wires ø 0.3 mm



CERAPALL® 6


White, copper-free, extra hard

Cerapall 6 has excellent thermal stability and is highly recommended for long span bridges. It is easy to cast and polish.
CTE 13.8 / 14.1 10⁻⁶K⁻¹



Recommended solders
before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information
200004790 25 g casting ingots
400003112 5 laser wires ø 0.3 mm



CERAPALL® 2


White, silver-free, extra hard

Cerapall 2 is known for its precise, trouble-free castings. It has uniquely fine and uniform grain structure which gives unmatched stability and accuracy in handling ultra thin margins.
Cerapall 2 has excellent mechanical properties and is easy to polish.
CTE 13.7 / 14.2 10⁻⁶K⁻¹



Recommended solders
before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information
200004763 25 g casting ingots
400003113 5 laser wires ø 0.3 mm





CERADELTA® 2


White, copper-free, extra hard

Ceradelta 2 is an excellent all-round alloy. It is easy to cast, has a light oxide colour and is easy to pre- or post-solder. It has a very high hardness.

CTE 14.5 / 14.9 10⁻⁶K⁻¹

Recommended solders
before bonding: VS 1 A
after bonding: Vacu PF

Ordering information
200004793 25 g casting ingots
400007057 5 laser wires ø 0.3 mm



CERADELTA®


White, copper-free, extra hard

Ceradelta is an alloy for long span bridges and partial cast dental frames with excellent properties for milling. It has good thermal stability and a light oxide colour.

CTE 15.1 / 15.3 10⁻⁶K⁻¹

Recommended solders
before bonding: VS 1 A
after bonding: Vacu PF

Ordering information
200004761 25 g casting ingots



METALLOY CC


White, extra hard

Metalloy CC has a low coefficient of thermal expansion. This ensures superb metal/ceramic bond behaviour and enables easy and extremely effective porcelain veneering. Metalloy CC also offers excellent laser welding ability since it is carbon free, and time is saved by the fact that oxide firing is not required.

CTE 14.0 10⁻⁶K⁻¹

Recommended solders
before bonding: Meta CC S1
after bonding: Meta CC S2

Ordering information
300005930 1000 g
300005931 500 g
300005932 100 g





AUROFLUID® PLUS

Yellow, palladium-free and copper-free, type IV, extra hard

Aurofluid Plus an outstanding high gold casting alloy with excellent casting, working and polishing properties. Especially designed for low fusing ceramics with high expansion, and has a low oxide formation.


CTE 15.8 / 16.2 10⁻⁶K⁻¹

Recommended solders

before bonding: Standard LFC 1 / Standard S 2

Ordering information

200004848 25 g casting ingots
400004405 5 laser wires ø 0.3 mm



PONTOR® LFC

Yellow, palladium-free, type IV, extra hard

Pontor LFC was developed as an universal alloy for conventional crowns and bridges with low fusing, high expansion ceramics.


CTE 16.7 / 17.0 10⁻⁶K⁻¹

Recommended solders

before bonding: Standard LFC 1 / Standard S 2
after bonding: Vacu 2

Ordering information

200004796 25 g casting ingots
400003145 5 laser wires ø 0.3 mm



PONTOR® 4CF

Light yellow, copper-free, type IV, extra hard

Pontor 4 CF is suitable for crowns, long span bridges, precision milling, partial dentures and low fusing ceramics. Pontor 4 CF is recommended for telescopic crowns.


CTE 16.2 / 16.8 10⁻⁶K⁻¹

Recommended solders

before bonding: Standard LFC 1 / Standard S 2
after bonding: Vacu 2

Ordering information

200004773 25 g casting ingots





AUROFLUID® 2PF

Yellow, palladium-free, type III, hard

Aurofluid 2 PF possesses a rich yellow colour. It has excellent finishing properties and polishes very easily. Aurofluid 2 PF is the classical inlay alloy and is also suitable for crowns and short span bridges.




Recommended solders

Vacu-PF / Standard S 2 / Standard S 4 PF

Ordering information

200004803 25 g casting ingots
400003154 5 laser wires ø 0.3 mm



PONTOR® MPF

Yellow, palladium-free, type IV, extra hard

Pontor MPF possesses a deep yellow colour. It is suitable for a variety of applications, including crowns, bridges, inlays and partial dentures. It has excellent mechanical properties.




Recommended solders

Vacu-PF / Standard S 3 / Standard S 4 PF

Ordering information

200004781 25 g casting ingots
400003103 5 laser wires ø 0.3 mm



AUROFLUID® 3

Yellow, type IV, extra hard

Aurofluid 3 is a conventional alloy with a yellow-reddish colour. It is suitable for inlays, crowns and bridges, milling work and cast partial denture-frameworks.




Recommended solders

Standard S 2 / Standard S 3 / Standard S 4 PF

Ordering information

200004733 25 g casting ingots
400003147 5 laser wires ø 0.3 mm





PONTOR® 2

Yellow, extra hard

Pontor 2 is an extra hard, medium gold content casting alloy. It handles beautifully and is recommended for crowns and long span bridges and precision milling works.




Recommended solders

Standard S 2 / Standard S 3 / Standard S 4 PF

Ordering information

200004758 25 g casting ingots
400003130 5 laser wires ø 0.3 mm



SOLARO® 3

Yellow, type IV, extra hard

Solaro 3 is an economical casting alloy with a yellow colour. It is suitable for long span bridges and milling works.




Recommended solders

Standard S 2 / Standard S 3 / Standard S 4 PF

Ordering information

200004736 25 g casting ingots
400003144 5 laser wires ø 0.3 mm



YELLOW SPECIAL

Light yellow, extra hard

Yellow Special is an economical alternative to high gold casting alloys. It provides excellent casting properties combined with ease of use and brings a yellow esthetic colour to the restorations.





Recommended solders

Standard S 3 / Standard S 4 PF / Standard S 31

Ordering information

200004844 25 g casting ingots





PAGALINOR® 2


White, hard

Pagalinator 2 is a white, hard casting alloy. It is suitable for inlays, crowns and short span bridges.



Recommended solders
Pagasol S 1 / Standard S 31

Ordering information
200004741 25 g casting ingots



PAGALIN® 2


White, hard

Pagalin 2 is an economical alloy due to its composition and its low density. It is suitable for inlays, crowns and short span bridges.



Recommended solders
Pagasol S 1 / Standard S 31

Ordering information
200004742 25 g casting ingots





IMPLANTALLOY 1


Yellow, extra hard

The rich golden colour of Implantalloy 1 assists in underlying all ceramic restorations with a correct hue. The light oxide layer enables excellent aesthetics even in situation with reduced space. Implantalloy 1 is recommended for all single crowns and three unit bridges in the anterior region.

CTE 13.8 / 14.1 10⁻⁶K⁻¹

Recommended solders
before bonding: VS 1 Plus
after bonding: Vacu PF

Ordering information
200016502 25 g casting ingots



IMPLANTALLOY 2


White, silver-free and copper-free, extra hard

The greatest emphasis was placed on the thermal stability for this silver-free and copper-free bonding alloy. Implantalloy 2 with its exceptionally high thermal stability shows excellent galvanic corrosion resistance behaviour when in contact with titanium.

CTE 13.8 / 14.1 10⁻⁶K⁻¹

Recommended solders
before bonding: VS 1 A / VS 1 C
after bonding: Vacu PF

Ordering information
200016503 25 g casting ingots




IMPLANTALLOY 3

Yellow, type IV, extra hard

With a Vickers hardness of 280 Implantalloy 3 is specially designed to offer ideal milling qualities. Despite its hardness Implantalloy 3 can be finished easily to the exacting standards required of all implant prosthetics.

Recommended solders
Standard S 2 / Standard S 3 / Standard S 4 PF

Ordering information
200016504 25 g casting ingots






CASTING




Precise casting of larger suprastructures is best achieved with indirect casting. We therefore recommend a transverse sprue of \varnothing 5mm, main feeder sprues of \varnothing 4mm and connectors of \varnothing 3mm; the connectors should be no longer than 4 mm.

Narrowing sprues in the contact area of the wax pattern encourages premature solidification in this area. This can lead to porosity within the casting, known as 'suck back' porosity.



Place the casting objects as far away as possible from the thermal centre of the ring and position them towards the outside to help directional freezing during casting.




The function of the casting reservoir is to supply the wax copings with liquid alloy during solidification. The casting reservoir should be placed as close to the thermal centre of the muffle as possible.

The quantity of alloy required to produce a perfect casting is calculated by weighing the waxed up copings, sprues and multiplying this by the specific gravity (density) of the alloy.



Preheating the casting crucible beforehand in the burnout furnace will reduce casting time and helps avoid overheating of the alloy.



The melting and casting of the alloy is one of the most critical stages. Take particular care to adjust the flame in the proper manner in order to achieve a perfect result. Keep the flame on the alloy constantly throughout the entire melting procedure. This provides additional protection for the alloy against oxidation.

Before commencing the melting procedure add a pinch of flux to the alloy. This addition of flux prevents oxidation of vital elements during melting.



SURFACE PREPARING



Prior to inspection of the fittings, remove any possible inclusions or flash under a binocular microscope.

The metal work is prepared in the usual manner. We recommend the use of tungsten carbide burs. Use only clean tungsten carbide burs. The use of 4 – 5 different tungsten carbide burs in shape and design is sufficient.




Finally, go over the metal work with a corundum stone at a maximum speed of 10,000 rpm.

Excessive pressure must be avoided at all times. Excessive pressure may damage the working parts of rotary instruments.

Use Cerasiv Brown grinding tips for high and medium gold bonding alloys and Cerasiv Blue for palladium based and non-precious bonding alloys.



Important: inspect the wall thickness of the metal work.



After final trimming of the metal work with a new corundum stone, oxidation may be carried out either directly, or the surface may be sandblasted prior to oxidation. Note that the copings must be cleaned by using ultrasound or a steam cleaner prior to oxidation.

The effect of sandblasting becomes visible after oxidation in the form of a more regular oxidation colour. The sandblasting procedure does not improve the adhesion of the opaque, as a rule.




SOLDERING



The area to be soldered should be flat and a minimum area of 6mm².

Use a rubber wheel to lightly polish the surface to be soldered.



The design of the solder block is a compromise between minimising its thermal mass whilst retaining sufficient strength to avoid its fracturing during soldering. Leave a parallel gap of 0.1 – 0.2 mm between surfaces to be soldered and sufficient area to ensure adequate strength of the joint.




Use a suitable instrument and apply the solder paste to the cold surfaces to be soldered as well as to the solder.

Heat the solder block to working temperature (approx. liquidus point of the solder).



Apply the solder at working temperature.

Should repeat applications of solder be necessary, only dip the solder into the solder paste. Do not apply more solder paste to the metalwork.





SOLDERING



Allow the solder block to bench cool.



Protect the ceramic with wax prior to contact with the soldering investment, for post bond soldering.



Keep the solder block as small as possible to ensure rapid and even heating. After steam cleaning, a flux is applied to the area to be soldered.



Dip the solder rods in the solder flux and crystallize them over the Bunsen burner. This minimises oxide formation on the solder and helps the soldering procedure.



After soldering in the furnace, remove the investment under running water. The oxide layer is removed using an etching agent. Trim the soldering areas with an abrasive rubber polishing wheel and give the restoration a final polish.



VS 1 C

White, silver-free and copper-free




Recommended solder for:

Implantalloy 2, V-Classic, V-Delta SF, V-Delta Special, V-Deltaloy
Cerapall 6, Cerapall 2

Working temperature: 1110°C
Melting range: 1030 - 1110°C


Ordering information

200004606 Ribbon on reel 3.0 g
400003134 Pack of 3 wires



VS 1 PLUS

Yellow, palladium-free and copper-free




Recommended solder for:

Implantalloy 1, V-Gnathos PF, V-Gnathos Plus, V-Gnathos Supra

Working temperature: 1040°C
Melting range: 980 - 1040°C


Ordering information

200004610 Ribbon on reel 3.0 g
400003131 Pack of 3 wires



VS 1 GPF

Yellow, palladium-free and copper-free




Recommended solder for:

V 92, V-Gnathos PF, V-Supragold

Working temperature: 1070°C
Melting range: 1010 - 1070°C

Ordering information

200004604 Ribbon on reel 3.0 g
400003133 Pack of 3 wires





VS 1 A

Light yellow, copper-free




Recommended solder for:

Ceradelta, Ceradelta 2, Cerapall 2, Cerapall 6, Implantalloy 2, V 92, V-Classic, V-Delta 450, V-Delta SF, V-Delta Special, V-Deltaloy

Working temperature: 1100°C
Melting range: 1030 - 1100°C

Ordering information

200004605 Ribbon on reel 3.0 g
400003132 Pack of 3 wires



VACU PF

Yellow, palladium-free




Recommended solder for:

All metalordental bonding alloys.

Working temperature: 850°C
Melting range: 750 - 800°C

Ordering information

200004727 Pack of rods 5.0 g



VACU 2

Light yellow, palladium-free




Recommended solder for:

Pontor 4 CF, Pontor LFC

Working temperature: 770°C
Melting range: 645 - 725°C

Ordering information

200004690 Pack of rods 5.0 g





STANDARD LFC 1

Yellow, palladium-free

Recommended solder for:

Aurofluid Plus, Pontor 4 CF, Pontor LFC




Working temperature: 880°C

Melting range: 855 - 880°C

Ordering information

200004608 Ribbon on reel 3.0 g

400003136 Pack of 3 wires





VACU-PF


Yellow, palladium-free

Recommended solder for:
Aurofluid 2 PF, Pontor MPF



Working temperature: 800°C
Melting range: 750 - 800°C


Ordering information
200004612 Ribbon on reel 3.0 g
400003207 Pack of 3 wires



STANDARD S 4 PF


Light yellow, palladium-free

Recommended solder for:
Aurofluid 2 PF, Aurofluid 3, Implantalloy 3, Pontor 2, Pontor 4 CF, Pontor MPF, Solaro 3, Yellow Special



Working temperature: 730°C
Melting range: 645 - 730°C

Ordering information
200004603 Ribbon on reel 3.0 g
400003139 Pack of 3 wires



STANDARD S 2


Yellow

Recommended solder for:
Aurofluid 2 PF, Aurofluid 3, Aurofluid Plus, Implantalloy 3, Pontor LFC, Pontor 2, Pontor 4 CF, Solaro 3



Working temperature: 830°C
Melting range: 790 - 830°C

Ordering information
200004601 Ribbon on reel 3.0 g
400003137 Pack of 3 wires





STANDARD S 3

Yellow




Recommended solder for:

Aurofluid 3, Implantalloy 3, Pontor 2, Pontor MPF, Solaro 3, Yellow Special

Working temperature: 780°C
Melting range: 710 - 780°C

Ordering information

200004602 Ribbon on reel 3.0 g
400003138 Pack of 3 wires



STANDARD S 31

Light yellow, palladium-free




Recommended solder for:

Pagalin 2, Pagalinor 2, Yellow Special

Working temperature: 790°C
Melting range: 775 - 790°C

Ordering information

200004613 Ribbon on reel 3.0 g



PAGASOL S 1

White




Recommended solder for:

Pagalin 2, Pagalinor 2

Working temperature: 850°C
Melting range: 800 - 850°C

Ordering information

200004607 Ribbon on reel 3.0 g
400003141 Pack of 3 wires





REVETOR II

Yellow, palladium-free

Recommended solder for:


Precious alloys to non-precious alloys.
May also be used for gilding clasps.



Working temperature: 800°C
Melting range: 750 - 800°C

Ordering information

- 200004592 Pack of 1 tube
- 200004593 Pack of 2 tubes



SOLCHROM PAL

White

Recommended solder for:


Precious alloys to non-precious alloys.



Working temperature: 940°C
Melting range: 865 - 940°C

Ordering information

- 200004590 Pack of 1 tube
- 200004591 Pack of 2 tubes






META CC S1


White

Recommended solder for:
Pre-solder for Metalloy CC



Working temperature: 1160°C


Ordering information
300005934 Pack of 3.0 g wire



META CC S2


Light yellow

Recommended solder for:
Post-solder for Metalloy CC



Working temperature: 720°C

Ordering information
300005935 Ribbon 0.4 g





ANCROFIL


Light yellow, palladium-free

For wires, clasps, attachments and construction elements.

Round wire, length = 200 mm

Ordering information

400003169	ø 0.40	1 x 0.40 g
400003172	ø 0.70	1 x 1.25 g
400003173	ø 0.80	1 x 1.65 g
400003175	ø 0.90	1 x 2.05 g
400003177	ø 1.00	1 x 2.55 g
400003179	ø 1.00	5 x 2.55 g
400003180	ø 1.10	1 x 3.05 g
400003182	ø 1.10	5 x 3.05 g
400003183	ø 1.20	1 x 3.65 g
400003185	ø 1.20	5 x 3.65 g
400003186	ø 1.30	1 x 4.25 g
400003170	ø 1.40	1 x 4.95 g
400003188	ø 1.50	1 x 5.70 g
400003189	ø 1.60	1 x 6.45 g
400003352	ø 1.80	1 x 8.15 g



NOVOSTIL

White, silver-free and copper-free

Non oxydizing dental alloy for the direct cast-on technique.


High melting range, 1400 - 1460°C.

For wires, posts, abutments, attachments and construction elements.

Round wire, length = 200 mm

Ordering information

400003190	ø 0.50	1 x 0.70 g
400003191	ø 0.70	1 x 1.40 g
400003192	ø 0.80	1 x 1.90 g
400003193	ø 0.90	1 x 2.30 g
400003194	ø 1.00	1 x 2.85 g
400003195	ø 1.10	1 x 3.45 g
400003196	ø 1.20	1 x 4.10 g
400003197	ø 1.30	1 x 4.80 g
400003198	ø 1.40	1 x 5.60 g
400003199	ø 1.50	1 x 6.40 g
400003201	ø 1.80	1 x 9.20 g
400003202	ø 2.00	1 x 11.40 g
400003203	ø 3.70	1 x 38.95 g





AUROFILM® 2000

Aurofilm 2000 is a preformulated gilding paste. Metalordental's expertise and unique experience in gold powder manufacturing has resulted in the development of a product with unique properties for the gilding of your crowns and bridges. It is easy to use, gives perfect covering on precious metals and a faultless appearance for your ceramic-to-precious metal restorations.

- Bright, natural looking shades on your metallo-ceramic restorations
- Golden base colour on all Pd- and Pd-Ag- Base alloys
- Easy to apply
- Ready to use
- Stable for multiple firings
- Excellent melting characteristics

Ordering information

200004585	Complete kit
200004586	3 g paste, syringe
200004393	Liquid 5 ml
200004492	Brush



AUROFILM® NP

Aurofilm NP has been developed for use on non-precious alloys as a gilding paste. Only one application is necessary prior to porcelain application.

Aurofilm NP is recommended to increase esthetics of Co-Cr and Ni-Cr non precious alloys.

Ordering information

200004588	Complete kit
200004589	3 g paste, syringe
200004394	Liquid 5 ml
200004492	Brush






HINTS FOR SUCCESS WITH THE METEOR FLAME-SYSTEM

Type of alloy	Type of burner	Heat soaking temp.	Type of crucible	Flux	Type of gaz
Conventional precious alloys	Meteor A	650°C - 700°C	preheat ceramic crucible in a furnace	Borax for alloys with a liquidus temperature below 1000°C Cristallized boric acid for alloys with a liquidus temperature above 1000°C	Propane / Compr. air
Universal precious alloys	Meteor O	650°C - 700°C	preheat ceramic crucible in a furnace	Borax for alloys with a liquidus temperature below 1000°C Cristallized boric acid for alloys with a liquidus temperature above 1000°C	Propane / Oxygen
Bonding precious alloys	Meteor O	750°C - 850°C	preheat ceramic crucible in a furnace	Cristallized boric acid	Propane / Oxygen
Non precious alloys	Meteor O-NP	900°C - 1000°C	preheat ceramic crucible in a furnace	none	Propane / Oxygen


Meteor A







Meteor O



Meteor O-NP



Pressure ratios	Ignition	Approximate length of flame	Shut down procedure	Neutral flame settings
Propane 0.5 - 1.0 bar Compr. air 1.0 - 2.0 bar	First open compr. air (blue), and then propane (red)	15 - 25 cm	Always close gas first (red), and then shut down compr. air (blue)	 Pointed blue cone up to mid length of a roaring flame
Propane 0.2 - 0.7 bar Oxygen 1.0 - 1.5 bar	First open oxygen (blue), and then propane (red)	20 - 22 cm	Always close gas first (red), and then shut down oxygen (blue)	 Soft drifting blue flame, inner flame cone 2 - 3 mm
Propane 0.2 - 0.7 bar Oxygen 1.0 - 1.5 bar	First open oxygen (blue), and then propane (red)	35 - 40 cm	Always close gas first (red), and then shut down oxygen (blue)	 Roaring hard blue flame, inner flame cone 3 - 4 mm
Propane 0.5 - 1.0 bar Oxygen 2.5 - 3.0 bar	First open oxygen (blue), and then propane (red)	40 - 45 cm	Always close gas first (red), and then shut down oxygen (blue)	 Hard roaring blue flame

Meteor L



Universal Handle





OXIDATION COLOUR OF METALDENTAL ALLOYS

Pontor® 4 CF, V 92, V-Gnathos® Plus, V-Gnathos® Supra, V-Supragold	
V-Gnathos® PF, Implantalloy 1	Cerapall® 2, V-Classic, V-Delta® 450
Aurofluid® Plus	V-Deltaloy
	Cerdelta® 2, Cerapall® 6, Implantalloy 2, Pontor® LFC, V-Delta® SF, V-Delta® Special
	Cerdelta®

Oxidation of the metal framework - and what this colour tells us:

The colour of the framework, after the recommended oxidation cycle, is a definite gauge as to the quality of the alloy surface. A regular oxide colour, specific to that alloy, indicates the optimal surface finish, clean and correctly processed, that provides a sound bond to the selected ceramic.

In order to achieve a successful reconstruction it is important not to underestimate the importance of the oxide layer as any discrepancies from the colour sample, described overleaf, points towards an incorrect treatment of the alloy.

AFFILIATES

Metalor Dental (Deutschland) GmbH

Augustenstrasse 124
D-70197 Stuttgart
Telefon +49 (0) 1803 638 2567
Fax +49 (0) 1803 638 2329
www.metalordental.de

Metalor Dental (France) SAS

Immeuble Le Sirius
9, rte du Colonel Marcel Moraine
F-92360 Meudon
Telefon +33 (0)1 46 99 95 60
Fax +33 (0)1 46 99 95 61
www.metalordental.com

Metalor Dental (Iberoamérica) S.A.

Avda. Diagonal 468, 5° A
E-08006 Barcelona
Telefon +34 933 032 311
Fax +34 933 074 707
www.metalordental.com

Metalor Dental Products (UK) Ltd.

Suite 10, The Green
Fountain Street, Macclesfield
Cheshire SK10 1JN
Telefon +44 (0)845 070 3595
Fax +44 (0)162 542 3558
www.metalordental.com

Metalor Dental (Hong Kong) Ltd.

Suite 1705-9, The Metropolis Tower
10 Metropolis Drive
Hung Hom
Kowloon, Hong Kong
Telefon +852 2521 4131
Fax +852 2365 1711
www.metalordental.com

Metalor Dental (Korea) Ltd.

8F. Shinhan Bldg. 2-42
Yangjae 1-dong
Seocho-gu
Seoul, Korea
Telefon +82 (0)2 575 3848
Fax +82 (0)2 575 2641
www.metalordental.com

Metalor Dental AG

Schweiz / Suisse / Svizzera
Bittertenstrasse 15
CH-4702 Oensingen
Telefon +41 (0) 800 292 292
Fax +41 (0) 62 388 69 70
www.metalordental.com

Metalor Dental (USA) inc.

255 John L. Dietsch Blvd
P.O. Box 255
02761 North Attleboro, MA
USA
Telefon +1 800 554 5504
Fax +1 508 695 3447
www.metalordental.com

HEADQUARTERS

Metalor Dental AG

Switzerland
Bittertenstrasse 15
CH-4702 Oensingen
Phone +41 (0)62 388 69 69
Fax +41 (0)62 388 69 70
www.metalordental.com

